


PRODUCT DATA SHEET

Revision 1	Pasteur pipette	
Validfor Item-No.:	3081000/3082000/3083000/3081001/3082001/3083001	

1.	Description / Specification	
1.1	Description	Pasteur pipette 3081000: Pasteur pipette, 1 ml, length 68 mm, bagged, non-sterile 3082000: Pasteur pipette, 2 ml, length 150 mm, bagged, non-sterile 3083000: Pasteur Pipette, 3 ml, 155 mm in length, bagged, non-sterile 3081001: Pasteur pipette, 1 ml, length 68 mm, bagged, sterile 3082001: Pasteur pipette, 2 ml, length 150 mm, bagged, sterile 3083001: Pasteur pipette, 3 ml, 155 mm in length, bagged, sterile
1.2	Dimensions	see Customer Drawing
1.3	Volume	3081000/3081001: 1ml 3082000/3082001: 2ml 3083000/3083001: 3ml
1.4	Material / Resin	Dish and Lid: PE (Polyethylene) The materials for manufacturing are <i>USP Class VI</i>
1.5	Colour	Dish and Lid: clear
1.6	Sterilization	3081000/3082000/3083000: non-sterile 3081001/3082001/3083001: sterile
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for bacteriology to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	100
3.2	Pieces / Box	3000
3.3	Lot-No.	DDD YY XXX (day, year, product No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
4.1	Laboratory use only

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 1 June 2024	Date 1 June 2024	Date 1 June 2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.

