


# PRODUCT DATA SHEET

Revision 1	200 µL Pipet tips	
Validfor Item-No.:	33011200/3011240/3011201/3011241/3011100/3011140/ 3011101/3011141	

1.	Description / Specification	
1.1	Description	200µL Pipet tips, 2~ 160µl  3011200 200µL Pipet tips, universal,yellow, non-sterile, bagged 3011240 200µL Pipet tips, universal,yellow, filter, non-sterile, bagged 3011201 200µL Pipet tips, universal,yellow, sterile, boxed 3011241 200µL Pipet tips, universal,yellow, filter, sterile, boxed 3011100 200µL Pipet tips, low retention, clear, non-sterile, bagged 3011140 200µL Pipet tips, low retention, clear, filter, non-sterile, bagged 3011101 200µL Pipet tips, low retention, clear, sterile, boxed 3011141 200µL Pipet tips, low retention, clear, filter, sterile, boxed
1.2	Dimensions	<u>Flask</u> : see Customer Drawing
1.3	Volume	Work volume: 2 ~ 160µl
1.4	Material / Resin	Pipette tip: PP (Polypropylene), free of heavy metal Filter: PE (Polyethylene) Rack/lid: PP (Polypropylene), free of heavy metal
1.5	Colour	Pipette tip: yellow/clear Filter: white
1.6	Sterilization	SAL 10 <sup>-6</sup>
1.7	Quality Control	<u>Raw Material-Control</u> : physical testing <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for pipetting of liquids to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only Expiry date and Lot-No. printed on the outer box

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -80°C to +121 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	3011200/3011100	3011240/3011140	3011201/3011241/ 3011101/3011141
3.1	Pieces / Bag	1000	1000	96
3.2	Pieces / Box	20	10	50
3.3	Lot-No.	DDD YY XXX (day, year, product No.)		
3.4	Other Information	Certificate of Quality to download		

4.	Other Information
4.1	Laboratory use only

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 1 June 2024	Date 1 June 2024	Date 1 June 2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.