


PRODUCT DATA SHEET

Revision 1	200µL Extended Pipet tips (77 mm)	
Validfor Item-No.:	3017000/3017010/3017001/3017011/3017100/3017110/ 3017101/3017111	

1.	Description / Specification	
1.1	Description	200µl Extended Pipet tips (77 mm), 2 ~ 300µl, clear 3017000 200µL Extended Pipet tips, universal,yellow, non-sterile, bagged 3017010 200µL Extended Pipet tips, universal,yellow, filter, non-sterile, bagged 3017001 200µL Extended Pipet tips, universal,yellow, sterile, boxed 3017011 200µL Extended Pipet tips, universal,yellow, filter, sterile, boxed 3017100 200µL Extended Pipet tips, low retention, clear, non-sterile, bagged 3017110 200µL Extended Pipet tips, low retention, clear, filter, non-sterile, bagged 3017101 200µL Extended Pipet tips, low retention, clear, sterile, boxed 3017111 200µL Extended Pipet tips, low retention, clear, filter, sterile, boxed
1.2	Dimensions	Flask: see Customer Drawing
1.3	Volume	Work volume: 2 ~ 200µl
1.4	Material / Resin	Pipette tip: PP (Polypropylene), free of heavy metal Filter: PE (Polyethylene) Rack/lid: PP (Polypropylene), free of heavy metal
1.5	Colour	Pipette tip: clear Filter: white
1.6	Sterilization	SAL 10 ⁻⁶
1.7	Quality Control	Raw Material-Control: physical testing Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for pipetting of liquids to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only Expiry date and Lot-No. printed on the outer box

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -80°C to +121 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	3017000/3017100	3017010/3017110	3017001/3017011/ 3017101/3017111
3.1	Pieces / Bag	1000	1000	96
3.2	Pieces / Box	20	10	50
3.3	Lot-No.	DDD YY XXX (day, year, product No.)		
3.4	Other Information	Certificate of Quality to download		

4.	Other Information
4.1	Laboratory use only

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 1 June 2024	Date 1 June 2024	Date 1 June 2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.